

Product Size	Color Code	Item #	Granule Size	Granule Net Volume	Content per Unit
Large		C11-008	500 - 1000µm	0.4 ml	3 systems
Medium		C11-078	500 - 1000µm	0.25 ml	3 systems
Small		C11-018	500 - 630µm	0.15 ml	3 systems




Available in three sizes, each unit contains three systems of the same size. A system of GUIDOR® *easy-graft*® CLASSIC Alloplastic Bone Grafting System contains one syringe of coated granules and one ampule of BioLinker®.

Size Selector – How much GUIDOR *easy-graft* is needed to fill a socket?

The following estimates are based on GUIDOR *easy-graft* placement in dentiform sockets. Material needs in clinical use may vary.

Maxillary



Mandibular